Module 3 | Lecture 1 Electrochemistry in Diagnostics

Instructor: Dr. Rajiv K. Kar

Jyoti and Bhupat Mehta School of Health Sciences and Technology

Indian Institute of Technology Guwahati




Biosensors | Point-of-Care Technology

- Refers to medical technologies that are designed for use near the patient, typically at the patient's bedside, rather than in a laboratory.
- Key Characteristics:
 Convenience: Portable, compact, and easy to use, enabling healthcare providers to perform tests quickly.
- * Rapid results: Provide results in minutes or hours.
- Patient comfort: Non-invasive or minimally invasive.
- Increased accessibility: Allows to perform tests and diagnose conditions in resource-limited settings | Cost effective solution.

Biosensors | Role of Electrodes


 A widely used method for sensing and analyzing | biological and chemical substances.

• An electrical signal is used to probe the presence of a target substance, and the resulting signal is measured and analyzed. 07/05/2023



3

Biosensors | Role of Electrodes



Source: https://doi.org/10.3390/s8042043


Electrochemical Analysis Potentiometric

Working Principle

- Uses the electrical potential difference between two electrodes to determine the concentration of a target substance.
- A reference electrode and a working electrode are immersed in a solution containing the target substance.
- A potential is applied across the electrodes, and the resulting potential difference is measured.
- The <u>potential difference</u> is proportional to the <u>log</u> of the concentration of the target substance, and the concentration can be determined from the measurement

Advantage:

- High accuracy, selectivity, and simplicity.
- Can be easily integrated into compact, portable ^{07/05/20}devices.




Electrochemical Analysis | **Amperometry**

Working Principle

- It measures the current produced by the oxidation or reduction of a target substance at an electrode.
- Consists of working electrode and a reference electrode, which are immersed in a solution containing the target substance.
- Applied potential causes the target substance to undergo oxidation or reduction at the working electrode.
- It generates a current that is proportional to the concentration of the target substance.
- The current can be measured and used to determine the concentration of the target substance.


Advantage:

High sensitivity, rapid response time, and the ability to detect trace amounts of target substances.



Amperometry vs Voltammetry

 Voltage Ramp: Voltammetry application of a linear voltage ramp to a working electrode, which is immersed in a sample solution. Leading to repeated oxidation and reduction cycles, generating a current that is proportional to the number of electroactive species.



Amperometry, on the other hand, is characterized by the application of a constant potential to the working electrode.

Current Measurement: Voltammetry measure the current generated as a function of time.

Amperometry measure the current generated by a sample at a single, constant potential.

 Information Obtained: Voltammetry -oxidation and reduction potentials of the electroactive species.


Amperometry, rate of electron transfer for a single, specific electrochemical reaction.


• Sample Characterization: Voltammetry is more useful for characterizing complex mixtures of electroactive species

Electrochemical AnalysisImpedimetric

Working Principle

- It measures the impedance (resistance to the flow of electrical current) of a system.
- A working electrode and a reference electrode, are immersed in a solution containing the target substance.
- Impedance of the system is determined by measuring the voltage and current that flow through the electrodes.
- Impedance of the system is influenced by several factors, including the presence of a target substance, the concentration of the o_{7/05/20}target substance, and the characteristics of the electrodes.





Z' (real; resistance) R_{ct}

high

frequency

WWW

low

frequency

Electrochemistry Fundamental

Working Principle

- Electrochemistry a tool to probe reactions involving electron transfers | flow of electrons to chemical changes.
- In inorganic chemistry, the resulting chemical change is often the oxidation or reduction of a metal complex.
- Example chemical reduction vs electrochemical reduction (reduction of ferrocenium [Fe(Cp)2] + (Cp = cyclopentadienyl), abbreviated as Fc+, to ferrocene [Fe(Cp)2], abbreviated as Fc)

Through a chemical reducing agent


$$Fc^+ + [Co(Cp^*)_2] \rightleftharpoons Fc + [Co(Cp^*)_2]^+$$


 Γ_{\bullet}^{\perp} , \bullet_{-} , Γ_{\bullet}

At an electrode

Source: http://dx.doi.org/10.1021/acs.jchemed.7b00


Electrochemistry Fundamental | **Graphics**






Cyclic Voltammetry

Oxidation High (V) Low Potentials



Working Principle

- The x-axis represents applied potential
 (E) and y-axis is the response i.e. curre
 (i) passed.
- The arrow indicates the beginning and sweep direction of the first segment (or ^{07/05/20} forward scan").



Source: http://dx.doi.org/10.1021/acs.jchemed.7b00

Next in the class

Cell Potential

SHE

Nernst Equation